

## Systems for Neutron Scattering

User Meeting May 2008




## Overview

- State of the art magnet technology for Neutron Scattering
- Low temperature sample environments for Neutron Scattering
- Remove or reduce dependency on supply of liquid helium

## The Design Process

- Finite Element Analysis
- Former Design
- Wire
- Quench Management
- Coil Structure




## The Design Process

- Finite Element Analysis
- Former Design
- Wire
- Quench Management
- Coil Structure

**PICTURE REMOVED**

## The Design Process

- Finite Element Analysis
- Former Design
- **Wire**
- Quench Management
- Coil Structure




## The Design Process

- Finite Element Analysis
- Former Design
- Wire
- Quench Management
- Coil Structure

**PICTURE REMOVED**

## The Design Process

- Finite Element Analysis
- Former Design
- Wire
- Quench Management
- Coil Structure



## Latest developments

- 6th 15T system supplied
  - Delivered within 6M of order (as per system 4 and 5)
- New horizontal field magnet design (6.5T field), with large 45 degree access
- Cryostat design improvements following consultation with users, in particular Eddy Lelievre of ILL
- 2 magnet systems to LLB, Saclay in 2008
  - 10T magnet for SANS
  - 10T asymmetric magnet (new technique for cadmium coating)
- Magnets under design/construction
  - 3 x unique designs, all with recondensing cryostats
  - 14T at 4K





## Overview

- State of the art magnet technology for Neutron Scattering
- Low temperature sample environments for Neutron Scattering
- Remove or reduce dependency on supply of liquid helium


## Pulse Tube Refrigerators

- PTR advantages over GM cryocooler designs
  - no moving parts in the low temperature region
  - significantly reduced vibration
  - The low temperature sections have no need for maintenance in normal use



## Recondensing Systems for Neutron Scattering

- ActivelyCooled™ magnet technology is a helium recondensing technology that reduces dependency on liquid cryogenics.
- Why recondensing and not Cryofree?
  - Stability of liquid bath for large complex coil structures
  - Improved vibration performance
  - System hold time during power failure
  - Integrated sample environments
  - Cool down before moving to beamline
  - Use conventional ULT inserts




## Cryofree™ sample environments – VeriCold VC4-Tn

- In close collaboration with the new neutron source FRM-2 in Garching, VeriCold Technologies GmbH has developed a pulse tube based 4K Closed Cycle Refrigerator for Neutron Scattering
  - Completely automatic
  - Initial cool down time of < 4 hours
  - Sample cool down time of < 2 hours from room temperature
  - Magnetic fields up to 9T can be applied
  - Windows for direct sample characterisation



## Cryofree™ sample environments – Heliox™ AC-V

- For wider temperature ranges down to 300 mK the **Heliox™ AC-V** offers Cryofree® operation at the touch a button.
  - Base temperatures < 300 mK
  - hold times > 50 h
  - High temperature operation up to 300 K
  - Thin wall Aluminum tails available
  - Seamless tails also available. Smooth high temperature
- This patented 3He technology developed on **Heliox™ AC-V** has led Oxford Instruments to develop our latest low temperature system **Triton™ DR**.





## Triton™ DR

- 50 (20) mK cryofree® base temperature, 10  $\mu$ W at 100 mK cooling power
- 300K operation
- High cooling power intermediate stage at 400 mK
- Low vibration pulse tube refrigerator technology
- No pumps or gas handling system
- Patented self-contained cryogenic cycle for leak-free reliable operation
- No needle valves or impedances




## Triton™ DR Operation



## VeriCold DR 200-10

- 10 mK cryofree® base temperature,
- 200  $\mu$ W at 100 mK cooling power
- JT stage provides cooling of the returning  $^3\text{He}$  mixture from 4K to below 1K.
- Oil free operation
- Single vacuum (No IVC, no indium seal)





## Summary

- For large complex superconducting magnets Actively Cooled™ recondensing technology offers stable magnet operation with minimum He consumption.
- Large split pair magnets are poor candidates for Cryofree® technology
- Low temperature sample environments are now available Cryofree®
- Triton DR Cryofree Dilution Refrigerators for Neutron Scattering